

COT Aging Report

All Experimenters Meeting

June 14, 2004

Morris Binkley and Kevin Burkett


Basics of Aging in COT

- Aging due to a coating on the sense wires.
- Analysis shows coating is mostly carbon and hydrogen with a little oxygen.
- Aging enhanced on exhaust side of COT.
- Presumably polymers (strings made of CH₂ groups?) grow in gas. When the strings get long enough, they plate out on the sense wires.
- Probably deposition rate is very temperature sensitive (temperature also controls convection).
- Probably are competing processes when avalanches are present: ablation and deposition.


Recirculation added to increase gas flow to limit polymer buildup.

- Started on May 14, 2004, but problems with the pump limited the flow.
- Full 160-170 SCFH established May 24, 2004.
- Since May 14th we have accumulated about 0.5 mC/cm on SL8 wires.
- Total Run2 accumulation on SL8 with HV on all Superlayers is about 6.9 mC/cm.
- Still not enough data to definitively measure aging with recirculation.
- Can monitor aging with currents and luminosity.
- Can monitor aging with widths of hits on tracks since widths are proportional to the charge in chamber pulses.


Cot Currents since April 2003




Width vs Date - SL1:4




Width vs Date - SL5:8



Width Relative To SL8 vs Date - SL1:4




Width Relative To SL8 vs Date - SL5:8




Can look COT currents on hourly basis



- Plots start after March shutdown and do not include first 6 hours of a store to reduce transient effects.
- First plot shows only currents for only those hours that inner super-layers are on:
 - Recirculation Starts May 14th (hour 42)
 - Full Recirculation May 24th (hour 90)
 - Cosmic Ray run May 27th (hour 119)
- Second plot only for outer super-layers and includes hours when inner super-layers are off (April and first part of May):
 - April 20th store with excess O₂ (hour 181-189).
 - Start of recirculation May 14 (hour 396) with excess O₂ from hour 384-417.



SL Current/Luminosity vs approx. hours HV on




SL Current Ratios vs approx hours HV on





Can check aging vs location in COT using widths of hits on good tracks.

- Aging depends of gas flow direction. Flow was reversed March 2, 2004.
- Inner Super-layers off all of April and half of May which allows chance to observe low aging.
- Divide COT in to six regions: two in z and three in ϕ and examine widths (which are proportional to charge in the chamber pulse).
- Recall that flow originally west to east.
- After reversal, flow is east to west.
- Pick upper (Φ_2) west as reference location.


Corrected Width SL2 vs. Date




SL2 Width Relative to Good Region vs. Date




Corrected Width SL6 vs. Date




Corrected Width SL6 vs. Date



SL6 Width Relative to Good Region vs. Date




SL6 Width Relative to Good Region vs. Date




Tests of bad wire plane using GMC with Sr90 Source


- Have tried with both Ar/Ethane/Alcohol and Ar/CO₂ gases.
- See reverse aging with both gases.
- See currents in chamber increasing linearly with time indicating there is a steady ablation of coating.
- Use Fe55 spectra to measure gain vs position.
- With Ar/Ethane/Alcohol the aging is reversed directly under the source but rest of chamber still ages.
- With Ar/CO₂ aging gets better everywhere.
- Following plot show effect of irradiating hole 7 of GMC. Black is before, red during, green after.




Fe55 Source Hole 8




Fe55 Source Hole 7




Fe55 Source Hole 6




Fe55 Source Hole 11




Fe55 Source Hole 9



Fe55 Source Hole 10



Fe55 Source Hole 13



Fe55 Source Hole 15

Conclusions

- Need more data, but aging probably improved with recirculation.
- Want to add O₂. Hopefully this week.
- If O₂ doesn't work will probably have two options:
 - 1) Try cleaning up the recirculated gas.
 - 2) Switch to Argon/CO₂ based gas
- Information at following web sites:
<http://ncdf82.fnal.gov/~binkley/cot/aging/>
<http://fcdfwww.fnal.gov/~burkett/COT/recent/>