T-1005 Report (g-2 Detector Team) April 23rd, 2012

Groups Involved (on site for all or part of the run)

Washington

Kentucky

Virginia


Fermilab

Shanghai

Cornell

Overview

- BNL E821 measured muon g-2 to 0.54 ppm
- Measurement was statistics limited
 - Increase statistics by moving to Fermilab
 - Goal > 21x increase
 - Reduce systematics with new detectors (among other things)
 - Goal ~ factor of 3 reduction
- Overall goal is factor of 4 improvement in measured value of muon g-2



Questions for this beam test

Calorimeter:

Tungsten SciFi or PbF2?

Photodetector:

Very fast PMT or SiPM?

Digitizer:

500 MHz or 1 GHz?

Equipment Deployed & Tested

EM Calorimeters:

- PbF2 array of 7 crystals in 2 x 3 x 2 stack (various wrapping, couplings)
- PbW04 crystal (as alternative in stack above)
- W/SciFi calorimeter
- PbGlass (Fermilab device as reference)

Readout devices

- Very fast Hamamatsu 8-Stage, R9800 PMT
- 16-Channel Hamamatsu SiPM 12 x 12 mm², 50 mm pitch, 56000 pixels, custom summed amplifier card
- Various Photonis and EMI 22 mm PMTs for side crystals

Beam defining equipment & trigger counters


- Fast MWPC with 3% Dimethoxy Methan 15% Isobutane, rest Argon
- 1.5 cm wide scintillator "finger counters"; larger paddle
- Cherenkov

Electronics / DAQ test equipment

- 8-channels of 12-bit, 500 MHz digitizers
- 2-channels of 8-bit 1 GHz digitizers
- MIDAS full system in test mode related to g-2 designed architecture


g-2 Program Goals and Highlights

- Short pulses
- Segmentation: 24 arrays of 5x7 channels
- Good resolution
- Dense (short radiation length)



PbF2 Cherenkov; Typically 3-4%/sqrt(E);

Pulse-optimized **SiPM array**, can be used in storage ring magnetic field Fast PMTs as backup; requires complex lightguides



W/SciFi; Typically 12-15%/sqrt(E). Very dense; slightly cheaper. Similar readout options




Setup




MIDAS DAQ

Crystal preparation

Fast PMT




SiPM




Preliminary Findings – PbF2

- Wrapping: Millipore paper; Black end; Optical grease
 - Variants: all white → longer pulse; all black → short pulse but low light; silicon interface → much lower light
- Pulse shape
 - FWHM ~8 ns; FWTM ~20 ns
- Resolution @ 4 GeV
 - 2.4% (includes $\Delta P/P$ of beam)



PbF2 Sum of 7 Crystals






Typical Program

- Calibrate all tubes at 4 GeV
- Sweep energies with open trigger
 - -1, 2, 3, 4, 5, 8 GeV
- Vary wrapping, couplings, detectors, readout devices, angles of impact, rate of beam
- Record pulse shapes, scope shots, long acquisition periods of digitized islands for pulse finding and energy reconstruction


16-channel SiPM — first tests Analysis of array in progress





12 x 12 mm²

W/SciFi – followup to our 2010 run (simple analysis)



The end ...



Thanks Fermilab for fantastic **SUPPORT** and excellent **BEAM** !!! *The g-2 Detector Team*