

Completion of ArgoNeuT NuMI Run

All Experimenter's Meeting


Mitch Soderberg

3/15/2010





ArgoNeuT

- ArgoNeuT is a ~175 liter Liquid Argon Time Projection Chamber (LArTPC) neutrino detector.
- Jointly funded by DOE/NSF
- Ran in NuMI beam at Fermilab, in front of MINOS near detector (to aid in muon reconstruction).
- Goals:
 - ▶ Gain experience building/running LArTPCs.
 - ▶ Accumulate neutrino/antineutrino events (1st time in the U.S., 1st time ever in a low-E beam).
 - ▶ Develop simulation/reconstruction for LArTPCs and compare with data.




Fermilab



MINOS Hall at Fermilab

ArgoNeuT: Collaboration



**6 Institutions,
20 collaborators**

F. Cavanna

University of L'Aquila

B. Baller, C. James, G. Rameika, B. Rebel
Fermi National Accelerator Laboratory

M. Antonello, R. Dimaggio, O. Palamara
Gran Sasso National Laboratory

C. Bromberg, D. Edmunds, P. Laurens, B. Page
Michigan State University

S. Kopp, K. Lang
The University of Texas at Austin


C. Anderson, B. Fleming, S. Linden, K. Partyka, M. Soderberg*, J. Spitz
Yale University

* = Spokesperson

LArTPC Principle: Reminder

TPC = Time Projection Chamber

- Neutrino interactions inside a TPC produce particles that ionize the argon as they travel (55k e⁻/cm).
- Ionization is drifted along E-field to wireplanes, consisting of wires spaced a few mm apart.
- Location of wires within a plane provides position measurements...multiple planes give independent views.
- Timing of wire pulse information is combined with drift speed to determine drift-direction coordinate.
- Scintillation light also present, can be collected by Photomultiplier Tubes and used in triggering.





Refs:

1.) *The Liquid-argon time projection chamber: a new concept for Neutrino Detector*, C. Rubbia, CERN-EP/77-08 (1977)

ArgoNeuT:TPC


Cryostat Volume	~500 Liters
TPC Volume	175 Liters
# Electronic Channels	480
Wire Pitch	4 mm
Plane Separation	4 mm
Electric Field	500V/cm
Max. Drift Time	330 μ s
Wire Properties	0.15mm diameter BeCu

Collection Induction #1 Induction #2



TPC About to Enter Cryostat (at PAB)

Used winding machines at Lab 6
Used routing machines at Lab 8




TPC Field Cage formed out of copper-clad G10 boards



Electric Field Rings

ArgoNeuT: Cryogenics

- Self-contained system....no refills.
- Continuously circulate argon through filters to purify.
- Cryocooler utilized to condense boil-off gas.
- Vacuum jacketed cryostat has 550 liter capacity.
- System designed to maintain ODH-0 rating of NuMI tunnel.




300W Cryocooler



ArgoNeuT: NuMI Run

- Filled the detector underground for the first time on Friday, May 8, 2009
- Acquired neutrino data for ~1 month before summer 2009 shutdown...continued running in the Fall, mostly in antineutrino mode
- Cryo. system operated continuously since initial fill, (modulo cryocooler repair for ~2 weeks in October).
- Run ended on Feb. 22, 2010, and system is now completely removed from the tunnel.





Moving underground
(lowering down 350 ft. shaft)



Installing underground.

ArgoNeuT: Removal



Safe and expedient process thanks to John Voirin and crew

ArgoNeuT: Removal




Safe and expedient process thanks to John Voirin and crew

ArgoNeuT: Removal




Safe and expedient process thanks to John Voirin and crew

ArgoNeuT Neutrino Event



Pixel size:
4mm x 0.3mm




ArgoNeuT Physics



- ArgoNeuT acquired $\sim 1.35\text{E}20$ Protons On Target (P.O.T.) by the end of its run, mostly in anti-neutrino mode.
- This data is being used to develop techniques for reconstructing events in 3D.
- Proving particle identification capability (dE/dx) using data will be an important result.
- We should obtain several cross-section measurements for the first time in a LAr experiment!

ArgoNeuT POT delivered and accumulated



Event Type	# in 180 days (1.4×10^{20} POT)
ν_{μ} CC	28800
$\bar{\nu}_{\mu}$ CC	2520
ν_e CC	540
NC	9720





Neutrino Mode

Event Type	# in 180 days (1.4×10^{20} POT)
ν_{μ} CC	9026
$\bar{\nu}_{\mu}$ CC	8111
ν_e CC	175
NC	5933

AntiNeutrino Mode

ArgoNeuT Software

- ArgoNeuT (anti)neutrino data spurring software/analysis work.
- We're creating our own LArTPC simulation/reconstruction software, nicknamed "LArSoft"
- LArSoft can/will be used for all future LAr experiments (MicroBooNE, LBNE, etc..).
- Example: Different reconstruction techniques being developed...



Conclusion


- ArgoNeuT run in the NuMI tunnel is now over...data analysis ramping up.
- Currently planning what to next with ArgoNeuT...stay tuned.

Thanks to everyone in Directorate, PPD, AD, ES&H, and CD who have helped make ArgoNeuT possible!

Back-Up Slides

Why Noble Liquids for Neutrinos?

- Abundant ionization electrons and scintillation light can both be used for detection.
- If liquids are highly purified (<0.1ppb), ionization can be drifted over long distances.
- Excellent dielectric properties accommodate very large voltages.
- Liquids are dense, so they make a good target for neutrinos.
- **Argon** is relatively cheap and easy to obtain (1% of atmosphere).




	He	Ne	Ar	Kr	Xe	Water
Boiling Point [K] @ 1 atm	4.2	27.1	87.3	120.0	165.0	373
Density [g/cm ³]	0.125	1.2	1.4	2.4	3.0	1
Radiation Length [cm]	755.2	24.0	14.0	4.9	2.8	36.1
dE/dx [MeV/cm]	0.24	1.4	2.1	3.0	3.8	1.9
Scintillation [γ /MeV]	19,000	30,000	40,000	25,000	42,000	
Scintillation λ [nm]	80	78	128	150	175	

Argon Purity

- Electron lifetime provides a measure of the purity of the liquid argon.
- We need a lifetime of several hundred microseconds in order to not lose an appreciable amount of charge as it drifts across the TPC.
- Plot below shows electron lifetime throughout the NuMI run.
- We learned a great deal about which cryogenic configurations gave the best purity improvements, and can use this information to achieve better performance in the future.

ArgoNeuT: Electron Lifetime and Impurities



ArgoNeuT: Electronics

- Bias voltage distribution & blocking capacitors on the TPC
- FET preamplifier similar to D0/ICARUS front-end
- Wide bandwidth filtering (10 - 159 kHz, now)
 - ▶ Full information on most hits/tracks
 - ▶ Employ DSP to extract hit/track parameters
- Digitization boards (ADF2 from D0) sample at 5 MHz (198ns), 2048 samples/channel
- Minimize noise sources
 - ▶ Double shielding of feed-through and preamplifiers
 - ▶ Remote ducted cooling
 - ▶ Extensive DC power filtering

